



2601B-PULSE Example TSP Script

DigitizeCompletePulse.tsp

Script Description

- This example script creates (and subsequently calls) a function that can be used with the Model 2601B-PULSE current pulser to output and digitize a single current pulse.
- With a few modifications, the basic approach shown in this example can be applied to the 2601B-PULSE SMU.
 - The integrating analog-to-digital converters must be used because the SMU does not have high-speed digitizers.
 - The pulser Pulse Width and Measure Delay attributes do not apply to the SMU. Instead Trigger Timers or SMU Source and Measure Delays must be used.

Note: Refer to the 2601B-PULSE Reference Manual for additional information about instrument operation and programming, including the Trigger Model.

Function Description

`digitize_pulse(current_level, pulse_width, meas_aperture, rangev, rangei, source_protectv, sense_protectv, bias_current)`

Pass Parameters:

- `current_level` : Peak current level of pulse in amps
- `pulse_width` : Width of current pulses in seconds
- `meas_aperture` : Effective integration time in seconds
- `rangev` : Voltage measure range in volts
- `rangei` : Current source and measure range in amps
- `source_protectv` : Voltage protection level at source terminals
- `sense_protect` : Voltage protection level at sense terminals
- `bias_current` : Idle current level in amps (base level for pulses)

Function Description - continued


- The pulse is generated using the instrument's Asynchronous Trigger Model.
- The Pulse Width is determined by a new pulser command, which provide more precise timing than a Trigger Timer. Needed to support pulses as short as 10 μ s.
 - `smua.trigger.source.pulsewidth`
- Using the dual 1MS/s digitizers built into the current pulser, the voltage and current are sampled simultaneously across the entire pulse.
 - Digitizers run at 1MS/s. Sample taken every 1 μ s.
 - Setting measure aperture (`smua.pulser.measure.aperture`) greater than 1 μ s causes multiple samples to be averaged and returned as a single reading.
 - Thus, the aperture setting defines an effective sample interval. If the aperture is set to 5 μ s, then a reading is taken every 5 μ s.
 - Therefore the effective sample rate (i.e. reading rate) is 1/measure aperture.

Function Description - continued


- The sampling nominally starts 10us before the start of the pulse and stops 10us after the end of the pulse.
- The pre-pulse acquisition interval is controlled using a standard Trigger Timer.
- The total acquisition time is the product of the measure count and the effective sample interval.
- See the next slide for a timing diagram.

Trigger Timing Scheme Used With Async Trig Model

Source Flow




Measure Flow



Function Description - continued

- Upon completion of the sweep, the data is printed to the Test Script Builder Instrument Console in a format that is suitable for copying and pasting into Microsoft Excel for graphing and analysis.
- Results are for: `digitize_pulse(10, 50e-6, 1e-6, 10, 10, 40, 20, 0)`



Time (s)	Current (A)	Voltage (V)
0.000000e+00	5.23664e-05	-4.54559e-05
1.000000e-06	-7.77895e-04	-5.29777e-04
2.000000e-06	-5.28838e-04	-3.68337e-04
3.000000e-06	-4.45809e-04	-4.54559e-05
4.000000e-06	-6.94876e-04	1.15985e-04
5.000000e-06	5.23664e-05	1.24607e-03
6.000000e-06	-2.79751e-04	9.23187e-04
7.000000e-06	-1.10997e-03	-4.54559e-05
8.000000e-06	-6.94876e-04	3.52643e-05
9.000000e-06	-7.77895e-04	-2.87617e-04
1.000000e-05	-4.45809e-04	-1.26176e-04
1.100000e-05	-4.45809e-04	3.58145e-04
1.200000e-05	2.70199e+00	2.91952e+00
1.300000e-05	8.16449e+00	9.06064e+00

Time (s)	Current (A)	Voltage (V)
0.000000E+00	5.23664E-05	-4.54559E-05
1.000000E-06	-7.77895E-04	-5.29777E-04
2.000000E-06	-5.28838E-04	-3.68337E-04
3.000000E-06	-4.45809E-04	-4.54559E-05
4.000000E-06	-6.94876E-04	1.15985E-04
5.000000E-06	5.23664E-05	1.24607E-03
6.000000E-06	-2.79751E-04	9.23187E-04
7.000000E-06	-1.10997E-03	-4.54559E-05
8.000000E-06	-6.94876E-04	3.52643E-05
9.000000E-06	-7.77895E-04	-2.87617E-04
1.000000E-05	-4.45809E-04	-1.26176E-04
1.100000E-05	-4.45809E-04	3.58145E-04
1.200000E-05	2.70199E+00	2.91952E+00
1.300000E-05	8.16449E+00	9.06064E+00

